

1262. CM + YE MEDIUM + 10% NaCl

Casamino acids Difco (vitamin free)	7.5	g
Yeast extract	10.0	g
MgSO ₄ x 7 H ₂ O	20.0	g
Na ₃ Citrat	3.0	g
KCl	2.0	g
NaCl	200.0	g
Fe*	10.0	ppm
Agar	15.0	g
NaCl (10%)	100.0	g
Distilled water	1000.0	ml

pH 7.4 (NaOH)

* Prepare 4.98% FeSO₄ x 7 H₂O (= 10 g Fe in 1l 0.001 M HCl); add 1 ml per liter medium.