

## 1279. HELICOBACTER MEDIUM

Brucella Broth (Becton Dickinson)	2.8	g
Aqua dem.	80.0	ml

Adjust with HCl to pH 5.0, make anoxic under a stream of N<sub>2</sub>/CO<sub>2</sub>, distribute in Hungate tubes to 8 ml and autoclave.

Add aseptically to each tube:

Fetal Bovine Serum (heat inactivated)	2.00	ml
Vitox (Oxoid)	0.06	ml
Skirrow supplement (Oxoid)	0.02	ml
Amphotericin (10 mg/ml)	0.02	ml

Add 1 ml of sterile air to each tube.