

1306. DEINOCOCCUS MEDIUM

Tryptone	5.0	g
Yeast extract	2.5	g
Dextrose	1.0	g
Maltose	1.0	g
DL-Methionine	1.0	g
Trace elements (see below)	1.00	ml
Vitamin solution (see below)	10.00	ml
Distilled water, refill to	1000.0	ml

Adjust pH to 7.0.

Vitamin solution:

Biotin	2.00	mg
Folic acid	2.00	mg
Pyridoxine-HCl	10.00	mg
Thiamine-HCl x 2 H ₂ O	5.00	mg
Riboflavin	5.00	mg
Nicotinic acid	5.00	mg
D-Ca-pantothenate	5.00	mg
Vitamin B ₁₂	0.10	mg
p-Aminobenzoic acid	5.00	mg
Lipoic acid	5.00	mg
Distilled water	1000.00	ml

Trace element solution SL-10:

HCl (25%; 7.7 M)	10.00	ml
FeCl ₂ x 4 H ₂ O	1.50	g
ZnCl ₂	70.00	mg
MnCl ₂ x 4 H ₂ O	100.00	mg
H ₃ BO ₃	6.00	mg
CoCl ₂ x 6 H ₂ O	190.00	mg
CuCl ₂ x 2 H ₂ O	2.00	mg
NiCl ₂ x 6 H ₂ O	24.00	mg
Na ₂ MoO ₄ x 2 H ₂ O	36.00	mg
Distilled water	990.00	ml

First dissolve FeCl₂ in the HCl, then dilute in water, add and dissolve the other salts. Finally make up to 1000.0 ml.