

1362. TRYPTONE SOYA BROTH (TSB) WITH ADDITIONAL SALTS

| | | |
|---------------------------------|--------|----|
| Peptone from casein | 17.0 | g |
| Peptone from soymeal | 3.0 | g |
| D(+)-Glucose | 2.5 | g |
| NaCl | 100.0 | g |
| K ₂ HPO ₄ | 2.5 | g |
| Distilled water | 1000.0 | ml |

Adjust pH to 7.2 – 7.5.

The basic medium is the same as Caso Bouillon (Merck or Oxoid). The medium may be solidified by adding 15.0 g/l agar.

When using commercially available Caso Bouillon (Merck or Oxoid) for strains DSM 24264, DSM 24313, DSM 24197 increase the NaCl concentration to 100.0 g/L.