

1366. GAUZE 1 WITH ADDITIONAL YEAST EXTRACT

Soluble Starch	20.00	g
NaCl	180.00	g
KNO ₃	1.00	g
K ₂ HPO ₄ x 3 H ₂ O	0.50	g
MgSO ₄ x 7 H ₂ O	0.50	g
FeSO ₄ x 7 H ₂ O	0.01	g
Yeast Extract	10.00	g
Distilled Water	1000.00	ml

Adjust to pH 7.0 and sterilize by autoclaving. The medium may be solidified by adding 20.0 g/l agar.