

1392. HALOTERRIGENA MEDIUM**Mineral Salts Solution:**

Casmino Acid	5.0	
NaCl	240.0	g
MgSO ₄ x 7 H ₂ O	10.0	g
KCl	5.0	g
Trisodium Citrate	3.0	g
KNO ₃	1.0	g
CaCl ₂ x 2 H ₂ O	0.2	g
Casmino Acid	5.0	g
Distilled Water	1000.0	ml

Adjust to pH 7.2 and sterilize by autoclaving. The medium may be solidified by adding 20.0 g/l agar.