

**1411. HALORUBRUM MEDIUM****Mineral Salts Solution:**

NaCl	195.00	g
MgCl <sub>2</sub>	32.50	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	50.80	g
CaCl <sub>2</sub> x 2H <sub>2</sub> O	0.80	g
KCl	5.00	g
NaHCO <sub>3</sub>	0.16	g
NaBr	0.60	g
Yeast extract	5.00	g
Distilled water	1000.00	ml

Adjust to pH 7.0 – 7.5 and sterilize by autoclaving. The medium may be solidified by adding 20.0 g/l agar.