

LENTIBACILLUS MEDIUM

NaCl	100.0	g
Casamino acids	5.0	g
Yeast extract	5.0	g
Glutamic acid	1.0	g
Trisodium citrate	3.0	g
MgSO ₄ x 7 H ₂ O	20.0	g
FeCl ₂ x 4 H ₂ O	36.0	mg
MnCl ₂ x 4 H ₂ O	0.4	mg
Agar, if required	20.0	g
Distilled water	1000.0	ml

Adjust to pH 7.2 and autoclave.