

1537. SPIRIBACTER MEDIUM (830 + 12% Marine salts)

R2A Broth	3.0	g
Marine Salts 12% (see below)	1000.0	ml
Agar	20.0	g
pH 7.0-7.2		

Prepare the medium. Sterilize by autoclaving. Before inoculation add from sterile stock solutions:

NaHCO ₃	0.71	g
--------------------	------	---

Marine Salts 12%:

NaCl	78.57	g
MgCl ₂ x 6 H ₂ O	15.37	g
MgSO ₄ x 7 H ₂ O	19.37	g
KCl	2.06	g
CaCl ₂ x 2H ₂ O	4.00	g
Distilled water	1000.00	ml