

1627. ACIDIBACTER MEDIUM

Glucose	0.90	g
Yeast extract	0.20	g
Salt solution	20.0	ml
Distilled water	980.00	ml

Adjust pH to 5.0 with sulfuric acid before autoclaving.

Salt solution

Na ₂ SO ₄	0.35	g
(NH ₄) ₂ SO ₄	2.25	g
KCl	0.25	g
MgSO ₄ x 7 H ₂ O	2.50	g
KH ₂ PO ₄	0.25	g
Ca(NO ₃) ₂ x 4 H ₂ O	0.07	g
Distilled water	100.00	ml