

1719. R2A WITH ARTIFICIAL SODA-SALT LAKE WATER

NaCl	46.94	g
Na ₂ SO ₄	7.84	g
MgCl ₂ x 6 H ₂ O	0.24	g
CaCl ₂ x 2 H ₂ O	0.17	g
KCl	1.38	g
KBr	0.19	g
H ₃ BO ₃	0.05	g
SrCl ₂ x 6 H ₂ O	0.08	g
NaF	0.06	g
Distilled water	900.0	ml

Add 3 g / L R2A (M. 830).

Autoclave to sterilize.

Add aseptically 100 mL sterile Na₂CO₃ solution (50 g / L, autoclaved).

pH 9.0.