

1755. A2 MEDIUM

| | | |
|--|----------|----|
| Sodium acetate | 0.240 | g |
| Yeast extract | 0.150 | g |
| MgSO ₄ x 7H ₂ O | 0.800 | g |
| CaCl ₂ x 2 H ₂ O | 0.030 | g |
| K ₂ HPO ₄ | 0.005 | g |
| KNO ₃ | 0.010 | g |
| HEPES | 20 mM | |
| Trace mineral salts solution | 2.000 | ml |
| a.dest | 1000.000 | ml |

Adjust pH to 7.0-7.5

Trace mineral salts solution

| | | |
|---|----------|----|
| FeCl ₃ x6H ₂ O | 1.000 | g |
| H ₃ BO ₃ | 2.500 | g |
| MnCl ₂ x4H ₂ O | 0.870 | g |
| CaCl ₂ x2H ₂ O | 3.700 | g |
| CuSO ₂ x5H ₂ O | 0.005 | g |
| ZnSO ₄ x7H ₂ O | 0.440 | g |
| Na ₂ MoO ₄ x2H ₂ O | 0.290 | g |
| ad distilled water | 1000.000 | ml |