

1796: PAS for *Acanthamoeba hatchetti*

Trisodium citrate x 2 H ₂ O	1.00	g
PAS Ah Stock Solution 1	5.00	ml
PAS Ah Stock Solution 2	5.00	ml
PAS Ah Stock Solution 3	5.00	ml

1. Make up to 1000 ml with MilliQ water.
2. Adjust pH to 7.0 with 1 N NaOH.
3. Autoclave at 121 °C for 20 min.
4. Store at room temperature.

PAS Ah Stock Solution 1

CaCl ₂ x 2 H ₂ O	5.90	g
--	------	---

Make up to 500 ml with MilliQ water.

PAS Ah Stock Solution 2

Na ₂ HPO ₄ (anhydrous)	35.50	g
KH ₂ PO ₄ (anhydrous)	34.10	g

Make up to 500 ml with MilliQ water.

PAS Ah Stock Solution 3

MgSO ₄ x 7 H ₂ O	98.60	g
--	-------	---

Make up to 500 ml with MilliQ water.