

1797: RE-103 Medium

NaCl, for Agar medium autoclave separately	206.00	g/l
MgSO ₄ x 7 H ₂ O	36.00	g/l
KCl	0.37	g/l
CaCl ₂ x 2 H ₂ O	0.50	g/l
Yeast extract (ThermoFischer BP1422)	5.00	g/l
MnCl ₂ x 2 H ₂ O	7.50	mg/l
Agar (optional)	15.00	g/l
Distilled water	1000.00	ml

1. Adjust pH to 7.0.
2. Note: When preparing a solid medium, autoclave the NaCl separately.