

## 1810: PSEUDOALTEROMONAS ESPEJIANA MEDIUM

Nutrient broth (Difco 234000)	8.00	g
KCl	0.70	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	1.50	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	12.00	g
NaCl	26.00	g
Agar (optional)	15.00	g
Distilled water	1000.00	ml

1. Adjust pH to 7.0.
2. Autoclave at 121 °C for 15 min.