

220. CASO AGAR (Merck 105458)

Peptone from casein	15.0	g
Peptone from soymeal	5.0	g
NaCl	5.0	g
Agar	15.0	g
Distilled water	1000.0	ml

Adjust pH to 7.3. Medium is identical with Tryptone Soya Agar (Oxoid Cm131).

For strain DSM 21449 und DSM 21450

NaCl	10.0	g/L
MnSO ₄	10.0	mg/L

For strain DSM 5746

Methanol	10.0	ml/L
----------	------	------