

223. POLYANGIUM MEDIUM

KNO ₃	1.00	g
K ₂ HPO ₄	1.00	g
MgSO ₄ x 7 H ₂ O	0.30	g
CaCl ₂ x 2 H ₂ O	0.13	g
FeCl ₃ x 6 H ₂ O	32.00	mg
Agar	10.00	g
Tap water	1000.00	ml

Adjust pH to 7.2 (use pH indicator-paper).

Aseptically add one or more sterile strips of filter paper to slants or plates. Inoculate the filter papers with a small (5 mm x 5mm) agar cube or a piece of filter paper overgrown with the strain.