

**512. MAGNETOSPIRILLUM II MEDIUM**

K <sub>2</sub> HPO <sub>4</sub>	0.5	g
Na-acetate	1.0	g
NH <sub>4</sub> Cl	0.1	g
Yeast extract	0.1	g
Fe(III) citrate	20.0	µM
Na-thioglycolate	0.5	g
Distilled water	1000.0	ml

Adjust pH to 6.8. Prepare medium 2 - 3 days before use. Dispense 5 ml of medium in sterile screw capped culture tubes. Inoculate with 10% (v/v) culture. Incubate tubes upright; do not shake.