

535b. REACTIVATION WITH LIQUID MEDIUM

Rehydrate and grow lyophilized cells from the ampoule in 5 ml liquid broth 545. Subsequent subculturing may be carried out in liquid medium 545 or on agar plates.

535. TRYPTONE SOYA AGAR (same as TRYPTICASE SOY BROTH AGAR)

545. TRYPTONE SOYA BROTH (TSB)

Peptone from casein	17.0	g
Peptone from soymeal	3.0	g
D(+)-Glucose	2.5	g
NaCl	5.0	g
K ₂ HPO ₄	2.5	g
Agar if required	15.0	g
Distilled water	1000.0	ml

Adjust pH to 7.3.

Medium 545 is the same as Caso Bouillon (Merck or Oxoid, Germany).