

545. TRYPTONE SOYA BROTH (TSB)

Peptone from casein	17.0	g
Peptone from soymeal	3.0	g
D(+)-Glucose	2.5	g
NaCl	5.0	g
K ₂ HPO ₄	2.5	g
Distilled water	1000.0	ml

Adjust pH to 7.3.

The same as Caso Bouillon (Merck or Oxoid, Germany).

For DSM 22112 dissolve ingredients, boil medium for 1 min., then cool to room temperature under 100% N₂ gas atmosphere. Dispense medium under same gas atmosphere in culture vessels and autoclave.