

738. BPHD Medium

NH ₄ Cl	1.00	g
Na ₂ HPO ₄ × 2 H ₂ O	0.42	g
Na ₂ H ₂ PO ₄ × H ₂ O	0.18	g
MgCl ₂ × 6 H ₂ O	0.10	g
CaCl ₂ × 2 H ₂ O	0.10	g
KCl	0.10	g
FeSO ₄ × 7 H ₂ O	0.04	g
Vitamin solution (see medium 141)	10.00	ml
Trace element solution (see medium 141)	10.00	ml
Cresol or phenol	1.00	mM*
Distilled water	1000.00	ml

Adjust pH to 6.5. Add cresol/phenol after autoclaving.

* (e.g.: 20 ml/l of a 4.7 g/l phenol-solution)