

798. TINDALLIA MEDIUM

Na ₂ CO ₃	8.0	g
KCl	0.2	g
NaCl	10.0	g
NH ₄ Cl	0.5	g
K ₂ HPO ₄	0.2	g
Yeast extract	4.0	g
Trace elements SL-10 (see medium 320)	1.0	ml
Vitamin solution (see medium 141)	1.0	ml
Na ₂ S x 9 H ₂ O	0.5	g
Distilled water	1000.0	ml

pH 9.0

Prepare the medium anaerobically without the vitamins. Boil the medium, cool under N₂ and dispense under N₂ and autoclave the medium. Add the vitamin solution to the cooled, autoclaved medium.