

812. SULFOBACILLUS DISULFIDOOXIDANS MEDIUM

$(\text{NH}_4)_2\text{SO}_4$	3.0	g
KCl	0.1	g
KH_2PO_4	0.5	g
$\text{MgSO}_4 \times 7 \text{ H}_2\text{O}$	0.5	g
$\text{Ca}(\text{NO}_3)_2 \times 4 \text{ H}_2\text{O}$	0.1	g
Yeast extract	0.1	g
Distilled water	1000.0	ml

Adjust pH to 2.25 prior to autoclaving. Before use add per liter 10 ml of a filter sterilized 10% (w/v) glutathione solution.