

954. MEDIUM FOR H. DOMBROWSKII

Casamino acids	5.00	g
Yeast extract (Oxoid)	5.00	g
TRIS	12.10	g
KCl	2.00	g
MgCl ₂ x 6 H ₂ O	20.00	g
CaCl ₂ x 2 H ₂ O	0.20	g
NaCl	200.00	g
Agar	20.00	g
Distilled water	1000.00	ml

Adjust pH to 7.4. Add the agar after dissolving all ingredients in the water and adjustment of pH.