

**956. MEDIUM FOR DSM 14457 AND 14458**

$\text{KH}_2\text{PO}_4$	2.000	g
$(\text{NH}_4)_2\text{SO}_4$	2.000	g
$\text{MgSO}_4 \times 7 \text{H}_2\text{O}$	0.125	g
NaCl	0.500	g
$\text{FeSO}_4 \times 7 \text{H}_2\text{O}$	0.002	g
Distilled water	1000.000	ml

pH 7.2. To the sterilised medium add sterile methanol (to a final concentration of 0.5%) or methylamine (to a final concentration of 0.3%).