

**276. HALOMONAS MEDIUM**

NaCl	80.00	g
Casamino acids (with vitamins)	7.50	g
Proteose peptone No. 3	5.00	g
Yeast extract	1.00	g
Na <sub>3</sub> -citrate	3.00	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	20.00	g
K <sub>2</sub> HPO <sub>4</sub>	0.50	g
Fe(NH <sub>4</sub> ) <sub>2</sub> (SO <sub>4</sub> ) <sub>2</sub> x 6 H <sub>2</sub> O	0.05	g
Agar, if necessary	15.00	g
Distilled water	1000.00	ml

Before autoclaving adjust pH to 7.5 with NaOH