

## 449. M17 MEDIUM FOR LACTIC STREPTOCOCCI

Peptone from casein	5.00	g
Soya peptone	5.00	g
Peptone bacteriological	5.00	g
Yeast extract	2.50	g
Ascorbic acid	0.50	g
MgSO <sub>4</sub> × 7 H <sub>2</sub> O	0.25	g
Na <sub>2</sub> -β-glycerophosphate	19.00	g
Lactose	5.00	g
Distilled water	1000.00	ml

Adjust pH to 6.9 ± 0.2.

Lactose should be sterilized by filtration.