

**465h. MINERALMEDIUM WITH 6-METHYLNICOTINATE**

Supplement medium 465 with selenite, tungstate, 2.5 - 5 g/l 6-methylnicotinate from a stock solution (autoclave separately) and 1 ml/l vitamin solution (see medium 461).

*Selenite/tungstate stock solution:*

NaOH	500.0	mg
Na <sub>2</sub> SeO <sub>3</sub> × 5 H <sub>2</sub> O	3.0	mg
Na <sub>2</sub> WO <sub>4</sub> × 2 H <sub>2</sub> O	4.0	mg
Distilled water	1000.0	ml

Autoclave separately. Add 1 ml per l medium.