

934. THERMOCOCCUS WAIOTAPUENSE MEDIUM

NaCl	1.500	g
Na ₂ SO ₄	4.000	g
KCl	0.700	g
NaHCO ₃	0.200	g
KBr	0.100	g
H ₃ BO ₃	0.030	g
MgCl ₂	5.800	g
CaCl ₂	1.500	g
SrCl ₂	0.025	g
Tryptone	5.000	g
Yeast extract	1.000	g
Na ₂ S x 9 H ₂ O	0.750	g
Resazurin	0.002	g
Sulphur	5 - 10.000	g
Distilled water	1000.000	ml

pH 7.0

Prepare the medium anaerobically under N₂. Steam the medium for 3 hours on three successive days.