

**1195. KUSHNERIA AURANTIA MEDIUM**

NaCl	78.00	g
MgCl <sub>2</sub> x 6 H <sub>2</sub> O	13.00	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	20.30	g
KCl	2.00	g
NaBr	0.23	g
NaHCO <sub>3</sub>	0.06	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	0.33	g
Yeast extract	5.00	g
Agar (for solid media)	20.00	g
Distilled water	1000.00	ml

Adjust to pH 7.2 – 7.4 and autoclave.