

## 1594. MEDIUM S FOR RHIZORHABDUS

Enzymatic casein hydrolysate	5.0	g
Glucose	2.5	g
.K <sub>2</sub> HPO <sub>4</sub> x 3 H <sub>2</sub> O	1.3	g
.KNO <sub>3</sub>	0.5	g
.MgSO <sub>4</sub> x 7 H <sub>2</sub> O	0.5	g
Agar, if necessary	13.0	g
Distilled water	1000.0	ml

Autoclave phosphate component separately and add after cooling. Adjust pH to 7.2.