

**1695. MARINES PLANCTOMYCES MEDIUM MODIFIED I**

Peptone (Bacto)	1.00	g
Yeast extract (Bacto)	1.00	g
HEPES	2.38	g
Mineral salt solution (see below)	20.00	ml
Trace element solution (see below)	1.00	ml
Artificial sea water (see below)	250.00	ml
Agar (Bacto)	15.00	g
<b>Double distilled water</b>	<b>730.00</b>	<b>ml</b>

Sterilize 20 min./121°C. After cooling to 60°C add:

N-Acetylglucosamine (10% stock solution)	10.0	ml
Glucose sol. (25%, sterile-filtered)	4.0	ml
Vitamin sol. (double conc.)(see below, medium 621)	5.0	ml

Adjust pH to 7.5 (the medium is only weakly buffered; one needs approx. 10 drops/l medium of 6 N KOH).

*Artificial sea water:*

NaCl	46.94	g
Na <sub>2</sub> SO <sub>4</sub>	7.84	g
MgCl <sub>2</sub> x 6 H <sub>2</sub> O	21.28	g
CaCl <sub>2</sub>	2.86	g
NaHCO <sub>3</sub>	0.38	g
KCl	1.38	g
KBr	0.192	g
H <sub>3</sub> BO <sub>3</sub>	0.052	g
SrCl <sub>2</sub>	0.08	g
NaF	0.006	g
<b>Double distilled water</b>	<b>1000.000</b>	<b>ml</b>

*Mineral salts solution:*

Nitrilotriacetic acid (NTA)	10.00	g
MgSO <sub>4</sub> x 7 H <sub>2</sub> O	29.70	g
CaCl <sub>2</sub> x 2 H <sub>2</sub> O	3.34	g
Na <sub>2</sub> MoO <sub>7</sub> O <sub>4</sub> x 2 H <sub>2</sub> O	12.67	mg
FeSO <sub>4</sub> x 7 H <sub>2</sub> O	99.00	mg
"Metals 44" (see below)	50.00	ml
Distilled water	950.00	ml

Dissolve the nitrilotriacetic acid, adjust the pH to 7.0 with KOH (about 7.3 g). Dissolve other salts separately, combine and adjust the pH to 6.8 with NaOH or H<sub>2</sub>SO<sub>4</sub>.

*"Metals 44":*

Na-EDTA	250.000	mg
ZnSO <sub>4</sub> x 7 H <sub>2</sub> O	1095.000	mg
FeSO <sub>4</sub> x 7 H <sub>2</sub> O	500.000	mg
MnSO <sub>4</sub> x H <sub>2</sub> O	154.000	mg
CuSO <sub>4</sub> x 5 H <sub>2</sub> O	39.200	mg
Co(NO <sub>3</sub> ) <sub>2</sub> x 6 H <sub>2</sub> O	24.800	mg
Na <sub>2</sub> B <sub>4</sub> O <sub>7</sub> x 10 H <sub>2</sub> O	17.700	mg
Distilled water	1000.000	ml

*Trace element solution:*

Na-Nitrilotriacetat	1.50	g
MnSO <sub>4</sub> x H <sub>2</sub> O	0.50	g
Fe SO <sub>4</sub> x 7 H <sub>2</sub> O	0.10	g
Co(NO <sub>3</sub> ) <sub>2</sub> x 6 H <sub>2</sub> O	0.10	g
ZnCl <sub>2</sub>	0.10	g
NiCl <sub>2</sub> x 6 H <sub>2</sub> O	0.05	g
H <sub>2</sub> SeO <sub>3</sub>	0.05	g
CuSO <sub>4</sub> x 5 H <sub>2</sub> O	0.01	g
AlK(SO <sub>4</sub> ) <sub>2</sub> x 12 H <sub>2</sub> O	0.01	g
H <sub>3</sub> BO <sub>3</sub>	0.01	g
Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	0.01	g
Na <sub>2</sub> WO <sub>4</sub> x 2 H <sub>2</sub> O	0.01	g
<b>Double distilled water</b>	1000.00	ml

*Vitamin solution (double conc.):*

Biotin	4.00	mg
Folic acid	4.00	mg
Pyridoxine-HCl	20.00	mg
Riboflavine	10.00	mg
Thiamine-HCl x 2 H <sub>2</sub> O	10.00	mg
Nicotinamide	10.00	mg
D-Ca-pantothenate	10.00	mg
Vitamin B <sub>12</sub>	0.20	mg
p-Aminobenzoic acid	10.00	mg
<b>Double distilled water</b>	1000.00	ml